


Refinements to the Mouse Vascular Access Button™

On October 22, 2019 we introduced an improved design of the magnetic one channel mouse Vascular Access Buttons and the protective aluminum caps. Part numbers and prices did not change. No changes to surgical procedures should be required.




BUTTON

Part nos: VABM1B/25, VABM1B/22, VABM1BSM/25, VABM1BSM/22



Old
.30 grams



New
.27 grams

smaller flange

gold plated magnet for durability


Base area, around which skin is closed, reduced by 40%

PROTECTIVE CAP

Part nos: VABM1C (red cap only)
VABM1C-5 (sold as package of 5—one of each color)



Old
.52 grams



VABM1C-5

New
.32 grams

New shape
40% lighter