

375/20PS

20ga single channel plastic swivel (pink), sterile



SPECIFICATIONS

Tubing size	20ga (.036x.024in, .91x.61mm)
Mates with	PE-60
Dead volume	12µL
Typical torque	.045 oz-in
Materials in fluid path	Stainless steel, Teflon
Dimensions (D x L)	.375x1.83in (9.5x46mm)
Weight	2.3g
Autoclavable	No
Compatible spring tethers	PS95, PS95H (.125in OD)
Typical application	Rat infusion, viscous compounds

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/swivels-mounts/swivels/plastic>

Instructions For Use https://www.instechlabs.com/hubfs/pdfs/manuals/375XXPS_IFU.pdf

Product Video

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com