

375/22PS

22ga single channel plastic swivel (blue), sterile



SPECIFICATIONS

Tubing size	22ga (.028x.016in, .71x.41mm)
Mates with	BTCOEX-22, VAHBPU-T22, PE-50
Dead volume	6µL
Typical torque	.035 oz-in
Materials in fluid path	Stainless steel, Teflon
Autoclavable	No
Dimensions (D x L)	.375x1.83in (9.5x46mm)
Weight	2.2g
Compatible spring tethers	PS95, PS95H (.125in OD)
Typical application	Rat infusion

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/swivels-mounts/swivels/plastic>

Instructions For Use https://www.instechlabs.com/hubfs/pdfs/manuals/375XXPS_IFU.pdf

Product Video

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com