

375/25PS

25ga single channel plastic swivel (purple), sterile



SPECIFICATIONS

Tubing size	25ga (.020x.010in, .51x.25mm)
Mates with	BTCOEX-25, VAHBPU-T25, BTPE-25
Dead volume	2.5µL
Typical torque	.020 oz-in
Materials in fluid path	Stainless steel, Teflon
Autoclavable	No
Dimensions (D x L)	.375x1.83in (9.5x46mm)
Weight	2.2g
Compatible spring tethers	PS62 (.085in OD)
Typical application	Mouse infusion

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/swivels-mounts/swivels/plastic>

Instructions For Use https://www.instechlabs.com/hubfs/pdfs/manuals/375XXPS_IFU.pdf

Product Video

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com