


Part Number: C10PU-MCA2422



not to scale

Full Description: Catheter for mouse carotid artery, soft 1Fr-2Fr, 8cm, sq tip. LS25, SP25/12.

Unit: each

Pouch Desc: Mouse CAC PU 1-2Fr

Pouch: 3x7.5

Approx Dead Vol: 9 μL

Connects to: 25ga

EO Sterilized: ☒

Intended Species: mouse

Vessel: carotid artery

Usage Notes:

TUBING

Tip: square

Length 1

L1: 2.5 cm

Material (IDxOD): PU 1Fr soft (.18x.35±.05mm)

Length 2

L2: 5.5 cm

Material (IDxOD): PU 2Fr (.43x.69±.03mm)

FEATURES

Distal (tip) end	Location	Fixed
Feature A:	a: cm	☐
Feature B:	b: cm	☐
Feature C: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 2.5 cm	☒

Proximal end		
Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1: LS25

Part 2: SP25/12