


Part Number: C10PU-MDD2404



not to scale

Full Description: Catheter for mouse duodenum, 1 to 2Fr, 31.5cm, square tip, collars @ 0.3, 0.5, 1.5, 2.5cm. Fits 25ga

Unit: each

Pouch Desc: Mouse DDC PU 1-2Fr 31.5cm

Pouch: 4x10

Approx Dead Vol: 9 μ L

Connects to: 25ga

EO Sterilized: ☒

Intended Species: mouse

Vessel: duodenum

Usage Notes:

TUBING

Tip: square

Length 1

L1: 1.5 cm

Material (IDxOD): PU 1Fr soft (.18x.35 $\pm$.05mm)

Length 2

L2: 30 cm

Material (IDxOD): PU 2Fr (.43x.69 $\pm$.03mm)

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 0.3 cm	☒
Feature B:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	b: 0.5 cm	☒
Feature C:	Long silicone collar for 2Fr (0.9mm OD), 3mm wide	c: 1.5 cm	☒
Proximal end			
Feature Z:		z: cm	☐
Feature Y:		y: cm	☐
Feature X:	Small silicone collar for 2Fr (0.9mm OD), 2mm wide	x: 29.0 cm	☒

PARTS ADDED TO POUCH

Part 1:

Part 2: