


Part Number: **C10PU-MFV1819**



not to scale

Full Description: Catheter for mouse femoral vein/artery, 1 to 3Fr, 10cm, collars @ 1.2, 5, 8.7, 9.0cm. Fits 22ga.

Unit: each **Pouch Desc:** Mouse FVC PU 1-3Fr 10cm **Pouch:** 3x7.5

Approx Dead Vol: 16 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: mouse **Vessel:** femoral vein or artery

Usage Notes:

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 5 cm **L2:** 5 cm
Material (IDxOD): PU 1Fr (.18x.35 \pm .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 \pm .03mm)

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 1.2 cm	<input type="checkbox"/>
Feature B:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	b: 5.0 cm	<input checked="" type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	z: 1.0 cm	<input checked="" type="checkbox"/>
Feature Y:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	y: 1.3 cm	<input checked="" type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**