


Part Number: C10PU-MJV2412



not to scale

Full Description: Catheter for mouse jugular vein, PU 1 to 3Fr, 15cm, square tip

Unit: each

Pouch Desc: Mouse JVC PU 1-3Fr 15cm

Pouch: 3x7.5

Approx Dead Vol: 30 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: mouse

Vessel: jugular vein

Usage Notes:

TUBING

Tip: square

Length 1

L1: 5 cm

Material (IDxOD): PU 1Fr soft (.18x.35 $\pm$.05mm)

Length 2

L2: 10 cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)

FEATURES

Distal (tip) end	Location	Fixed
Feature A:	a: cm	☐
Feature B:	b: cm	☐
Feature C: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 5.0 cm	☒

Proximal end		
Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: