


Part Number: **C10PU-RBD1423**



not to scale

Full Description: Catheter for rat bile duct, 1Fr to VAHBPU-T22, 45cm, collar @ 9mm. Fits 22ga.

Unit: each **Pouch Desc:** Rat BDC 1Fr-024x055 PU 45cm **Pouch:** 4x10

Approx Dead Vol: 130 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** bile duct

Usage Notes:

TUBING

Tip: square **Length 1** **Length 2**
L1: 5 cm **L2:** 40 cm
Material (IDxOD): PU 1Fr (.18x0.35 \pm .05mm) **Material (IDxOD):** VAHBPU-T22 (.61x1.40 \pm .05mm)

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 0.9 cm	<input type="checkbox"/>
Feature B:		b: cm	<input type="checkbox"/>
Feature C:	Extra large silicone collar (2.2mm OD), 3mm wide	c: 5.0 cm	<input checked="" type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**