


Part Number: **C10PU-RFA1702**



*not to scale*

**Full Description:** Catheter for rat femoral artery, 1 to 3Fr, 27.5cm, collar @2cm. Fits 22ga.

**Unit:** each      **Pouch Desc:** Rat FAC 1-3Fr 30cm      **Pouch:** 4x10

**Approx Dead Vol:** 75  $\mu$ L      **Connects to:** 22ga      **EO Sterilized:**

**Intended Species:** rat      **Vessel:** femoral artery

**Usage Notes:**

## TUBING

**Tip:** rounded      *Length 1*      *Length 2*  
**L1:** 2.5 cm      **L2:** 25 cm  
**Material (IDxOD):** PU 1Fr (.18x.35 $\pm$ .05mm)      **Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm)

## FEATURES

*Distal (tip) end*      *Location*      *Fixed*

**Feature A:** Small silicone collar for 1Fr (0.6mm OD), 1mm wide      **a:** 2.0 cm     

**Feature B:**      **b:**   cm     

**Feature C:** Large silicone collar for 3Fr (1.7mm OD), 3mm wide      **c:** 2.5 cm     

*Proximal end*

**Feature Z:**      **z:**   cm     

**Feature Y:**      **y:**   cm     

**Feature X:**      **x:**   cm     

## PARTS ADDED TO POUCH

**Part 1:**      **Part 2:**