


Part Number: **C10PUS-RBD1810**



not to scale

Full Description: Catheter for rat bile duct, soft 1Fr to 3Fr, 58cm, collars @ 0.4, 3.0cm. Fits 22ga.

Unit: each **Pouch Desc:** Rat BDC PU soft 1Fr 54cm **Pouch:** 4x10

Approx Dead Vol: 175 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** bile duct

Usage Notes:

TUBING

Tip: square **Length 1** **Length 2**
L1: 4 cm **L2:** 54 cm
Material (IDxOD): PU 1Fr soft (.18x.35 \pm .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 \pm .03mm)

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 0.4 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	b: 3.0 cm	<input type="checkbox"/>
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 4.0 cm	<input checked="" type="checkbox"/>

<i>Proximal end</i>			
Feature Z:		z: <input type="text"/> cm	<input type="checkbox"/>
Feature Y:		y: <input type="text"/> cm	<input type="checkbox"/>
Feature X:		x: <input type="text"/> cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**