


Part Number: C10SS-MTV2417P



not to scale

Full Description: Mouse tail vein cannula, 29ga, 65cm PE-10 tube, PinPort and injector

Unit: each

Pouch Desc: Mouse TVC w/PP 65cm

Pouch: 4x10

Approx Dead Vol: 40 μ L

Connects to:

EO Sterilized: ☒

Intended Species: mouse

Vessel: tail vein

Usage Notes:

TUBING

Tip: sharp

Length 1

Length 2

L1: 65 cm

L2:

Material (IDxOD): PE-10 (0.3x0.6mm)

Material (IDxOD):

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	29ga stainless steel needle (.18x.33mm)	a: 1.2 cm	☒
Feature B:		b: cm	☐
Feature C:		c: cm	☐

Proximal end			
Feature Z:		z: cm	☐
Feature Y:	PinPort, white, 28ga connector	y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1: PNP3M injector

Part 2: needle sleeve