


Part Number: **C15SS-RTV1453P**



not to scale

Full Description: Rat tail vein cannula, 25ga, 6cm PE/PVC tube, PinPort and injector

Unit: each **Pouch Desc:** Rat TVC w/PP 6cm **Pouch:** 3x7.5

Approx Dead Vol: 12 μ L **Connects to:** **EO Sterilized:**

Intended Species: rat **Vessel:** tail vein

Usage Notes:

TUBING

Tip: sharp *Length 1* *Length 2*
 L1: 6 cm L2: cm
 Material (IDxOD): BTCOEX-25 (.4x1.3mm) Material (IDxOD):

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	25ga stainless steel needle (.25x.51mm)	a: 2.0 cm	<input checked="" type="checkbox"/>
Feature B:		b: cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:	PinPort, white, 25ga connector	y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: PNP3M injector **Part 2:** needle guard