


Part Number: C19PU-GCA2305



not to scale

Full Description: Catheter for GP carotid artery, PU 1.9-3Fr 15cm, collar @ 2.0cm, square tip. Fits 22ga.

Unit: each

Pouch Desc: GP CAC 1.9-3Fr 15cm

Pouch: 3x7.5

Approx Dead Vol: 40 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: guinea pig

Vessel: carotid artery

Usage Notes:

TUBING

Tip: square

Length 1

L1: 2 cm

Material (IDxOD): PU 1.9Fr (.31x.64 $\pm$.05mm)

Length 2

L2: 13 cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)

FEATURES

Distal (tip) end	Location	Fixed
Feature A:	a: cm	☐
Feature B:	b: cm	☐
Feature C: Large silicone collar for 3Fr (1.7mm OD), 3mm wide	c: 2.0 cm	☒

Proximal end		
Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: