


Part Number: C19PU-GCA2409



not to scale

Full Description:
Catheter for GP carotid artery, PU 1.9-3Fr 17cm, collar @ 4.0cm, square tip. Fits 22ga.

Unit: each

Pouch Desc: GP CAC 1.9-3Fr 17cm

Pouch: 3x7.5

Approx Dead Vol: 42  $\mu$ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: guinea pig

Vessel: carotid artery

Usage Notes:

TUBING

Tip: square

Length 1
L1: 4 cm

Length 2
L2: 13 cm

Material (IDxOD): PU 1.9Fr (.31x.64 $\pm$ .05mm)

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$ .03mm)

FEATURES

Distal (tip) end

Feature A:

Feature B:

Feature C: Large silicone collar for 3Fr (1.7mm OD), 3mm wide

Location

a:  cm

b:  cm

c: 4.0 cm

Fixed

☐

☐

☒

Proximal end

Feature Z:

Feature Y:

Feature X:

z:  cm

y:  cm

x:  cm

☐

☐

☐

PARTS ADDED TO POUCH

Part 1:

Part 2: