


Part Number: **C19PU-RCA1301**



*not to scale*

**Full Description:** Catheter for rat carotid artery, PU 1.9-3Fr 16.5cm, collars @ 3.5, 8cm. Fits 22ga.

**Unit:** each **Pouch Desc:** Rat CAC 1.9-3Fr PU 16.5cm **Pouch:** 3x7.5

**Approx Dead Vol:** 40  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** carotid artery

**Usage Notes:**

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 3.5 cm **L2:** 13 cm  
**Material (IDxOD):** PU 1.9Fr (.31x.64 $\pm$ .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm)

## FEATURES

Distal (tip) end		Location	Fixed
Feature A:		a: <input type="text"/> cm	<input type="checkbox"/>
Feature B:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	b: 3.5 cm	<input type="checkbox"/>
Feature C:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	c: 8.0 cm	<input type="checkbox"/>
Proximal end			
Feature Z:		z: <input type="text"/> cm	<input type="checkbox"/>
Feature Y:		y: <input type="text"/> cm	<input type="checkbox"/>
Feature X:		x: <input type="text"/> cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**