


Part Number: **C19PU-RFV1414**



not to scale

Full Description: Catheter for rat femoral vein, PU 1.9Fr 4.5cm, 3Fr 25cm, collars @ 4.0, 4.5, 7.5cm. Fits 22ga.

Unit: each **Pouch Desc:** Rat FVC PU 1.9-3Fr 30cm **Pouch:** 4x10

Approx Dead Vol: 80 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** femoral vein

Usage Notes:

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 4.5 cm **L2:** 25 cm
Material (IDxOD): PU 1.9Fr (.31x.64 \pm .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 \pm .03mm)

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 4.0 cm	<input type="checkbox"/>
Feature B:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	b: 4.5 cm	<input type="checkbox"/>
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 7.5 cm	<input type="checkbox"/>
<i>Proximal end</i>			
Feature Z:		z: <input type="text"/> cm	<input type="checkbox"/>
Feature Y:		y: <input type="text"/> cm	<input type="checkbox"/>
Feature X:		x: <input type="text"/> cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**