


Part Number: **C19PU-RFV1719**



*not to scale*

**Full Description:** Catheter for rat femoral vein, PU 1.9Fr 4.5cm, 3Fr 16cm, collars @ 4.0, 4.5, 7.5cm. Fits 22ga.

**Unit:** each **Pouch Desc:** Rat FVC PU 1.9-3Fr 30cm **Pouch:** 4x10

**Approx Dead Vol:** 55  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** femoral vein

**Usage Notes:**

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 4.5 cm **L2:** 16 cm  
**Material (IDxOD):** PU 1.9Fr (.31x.64 $\pm$ .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm)

## FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
<b>Feature A:</b>	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 4.0 cm	<input type="checkbox"/>
<b>Feature B:</b>	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	b: 4.5 cm	<input type="checkbox"/>
<b>Feature C:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 7.5 cm	<input type="checkbox"/>

<i>Proximal end</i>			
<b>Feature Z:</b>		z: <input type="text"/> cm	<input type="checkbox"/>
<b>Feature Y:</b>		y: <input type="text"/> cm	<input type="checkbox"/>
<b>Feature X:</b>		x: <input type="text"/> cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**