


Part Number: **C19PU-RJV1436**



not to scale

Full Description: Catheter for rat jugular vein, PU 1.9-3Fr 16.5cm, collars @ 2.8, 3.5, 8cm. Fits 22ga.

Unit: each **Pouch Desc:** Rat JVC 1.9-3Fr PU 16.5cm **Pouch:** 3x7.5

Approx Dead Vol: 45 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** jugular vein

Usage Notes:

TUBING

Tip: rounded **Length 1** L1: 3.5 cm **Length 2** L2: 13 cm
Material (IDxOD): PU 1.9Fr (.31x.64 \pm .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 \pm .03mm)

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 2.8 cm	<input type="checkbox"/>
Feature B:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	b: 3.5 cm	<input checked="" type="checkbox"/>
Feature C:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	c: 8.0 cm	<input type="checkbox"/>

Proximal end		Location	Fixed
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**