


Part Number: **C20PU-MDD1812**



not to scale

Full Description: Catheter for mouse duodenum, PU 2Fr 10cm, beveled tip, collars @ 0.8, 1.0cm. Fits 25ga.

Unit: each **Pouch Desc:** Mouse DDC 2Fr PU 10cm **Pouch:** 3x7.5

Approx Dead Vol: 15 μ L **Connects to:** 25ga **EO Sterilized:**

Intended Species: mouse **Vessel:** duodenum

Usage Notes:

TUBING

Tip: beveled *Length 1* *Length 2*
L1: 10 cm **L2:** cm
Material (IDxOD): PU 2Fr (.43x.69 \pm .03mm) **Material (IDxOD):**

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 0.8 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	b: 1.0 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**