


Part Number: C20PU-MDD2403



not to scale

Full Description: Catheter for mouse duodenum, PU 2Fr 30cm, collars @ 0.3, 0.5, 1.2cm, square tip. Fits 25ga.

Unit: each

Pouch Desc: Mouse DDC 2Fr PU 30cm sq

Pouch: 4x10

Approx Dead Vol: 45 μ L

Connects to: 25ga

EO Sterilized: ☒

Intended Species: mouse

Vessel: duodenum

Usage Notes:

TUBING

Tip: square

Length 1

L1: 30 cm

Material (IDxOD): PU 2Fr (.43x.69 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end

	Location	Fixed
Feature A: Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 0.3 cm	☒
Feature B: Small silicone collar for 2Fr (0.9mm OD), 1mm wide	b: 0.5 cm	☒
Feature C: Small silicone collar for 2Fr (0.9mm OD), 1mm wide	c: 1.2 cm	☒

Proximal end

Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: