


Part Number: C20PU-MJV1410

*not to scale*

Full Description: Catheter for mouse jugular vein, PU 2Fr 6cm, collar @ 1.3cm. Fits 25ga.

Unit: each Pouch Desc: Mouse JVC 2Fr PU 6cm Pouch: 3x7.5

Approx Dead Vol: 9 μ L Connects to: 25ga EO Sterilized: ☒

Intended Species: mouse Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 6 cm

Length 2

L2: cm

Material (IDxOD): PU 2Fr (.43x0.69 $\pm$.03mm)

Material (IDxOD):

FEATURES

Distal (tip) end

Feature A: Small silicone collar for 2Fr (0.9mm OD), 1mm wide

Location

Fixed

a: 1.3 cm

☐

Feature B:

b: cm

☐

Feature C:

c: cm

☐

Proximal end

Feature Z:

z: cm

☐

Feature Y:

y: cm

☐

Feature X:

x: cm

☐

PARTS ADDED TO POUCH

Part 1:

Part 2: