


Part Number: C20PU-MJV1425

*not to scale*

Full Description: Catheter for mouse jugular vein, PU 2Fr 8cm, collar @ 1.0cm. Fits 25ga.

Unit: each Pouch Desc: Mouse JVC 2Fr PU 8cm Pouch: 3x7.5

Approx Dead Vol: 12 μ L Connects to: 25ga EO Sterilized: ☒

Intended Species: mouse Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded Length 1 L1: 8 cm Length 2 L2: cm

Material (IDxOD): PU 2Fr (.43x0.69 $\pm$.03mm) Material (IDxOD):

FEATURES*Distal (tip) end*

	Location	Fixed
Feature A: Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 1.0 cm	☐
Feature B:	b: cm	☐
Feature C:	c: cm	☐

Proximal end

Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1: Part 2: