


Part Number: C20PU-MJV1511

*not to scale*

Full Description: Catheter for mouse jugular vein, PU 2Fr 10cm, collars @ 1.3, 5.8 cm. Fits 25ga.

Unit: each Pouch Desc: Mouse JVC 2Fr PU 10cm Pouch: 3x7.5

Approx Dead Vol: 15 μ L Connects to: 25ga EO Sterilized: ☒

Intended Species: mouse Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 10 cm

Length 2

L2: cm

Material (IDxOD): PU 2Fr (.43x.69 $\pm$.03mm)

Material (IDxOD):

FEATURES

Distal (tip) end

Feature A: Small silicone collar for 2Fr (0.9mm OD), 1mm wide

Location

a: 1.3 cm

Fixed

☐

Feature B: Small silicone collar for 2Fr (0.9mm OD), 1mm wide

b: 5.8 cm

☐

Feature C:

c: cm

☐

Proximal end

Feature Z:

z: cm

☐

Feature Y:

y: cm

☐

Feature X:

x: cm

☐

PARTS ADDED TO POUCH

Part 1:

Part 2: