


Part Number: C20PU-MJV1715VAB

*not to scale***Full Description:** Catheter for mouse jugular vein, PU 2Fr 8cm, collar @ 1.1cm, sq tip, VABM1B/25 attached, injector**Unit:** each **Pouch Desc:** Mouse JVC 2Fr PU 8cm VAB **Pouch:** 3x7.5**Approx Dead Vol:** 15 μ L **Connects to:** 25ga **EO Sterilized:** ☒**Intended Species:** mouse **Vessel:** jugular vein**Usage Notes:****TUBING****Tip:** square*Length 1***L1:** 8 cm*Length 2***L2:** cm**Material (IDxOD):** PU 2Fr (.43x.69 $\pm$.03mm)**Material (IDxOD):****FEATURES***Distal (tip) end***Feature A:** Small silicone collar for 2Fr (0.9mm OD), 1mm wide*Location**Fixed***a:** 1.1 cm☐**Feature B:****b:** cm☐**Feature C:****c:** cm☐*Proximal end***Feature Z:** Mouse VAB, 1ch magnetic, 25ga, bonded to catheter**z:** cm☐**Feature Y:****y:** cm☐**Feature X:****x:** cm☐**PARTS ADDED TO POUCH****Part 1:** PNP3M injector**Part 2:**