


Part Number: C20PU-MJV1814VAB



Full Description: Catheter for mouse jugular vein, PU 2Fr 8cm, collar @ 1.2cm, round tip, VABM1B/25 attached, injector

Unit: each Pouch Desc: Mouse JVC 2Fr PU 8cm VAB Pouch: 3x7.5

Approx Dead Vol: 15 μ L Connects to: 25ga EO Sterilized: ☒

Intended Species: mouse Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1: L1: 8 cm

Length 2: L2: cm

Material (IDxOD): PU 2Fr (.43x.69 $\pm$.03mm) Material (IDxOD):

FEATURES

Distal (tip) end

Feature	Description	Location	Fixed
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 1.2 cm	☐
Feature B:		b: cm	☐
Feature C:		c: cm	☐

Proximal end

Feature	Description	Location	Fixed
Feature Z:	Mouse VAB, 1ch magnetic, 25ga, bonded to catheter	z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1: PNP3M injector Part 2: