


Part Number: C20PU-MJV1818VAB



Full Description: Catheter for mouse jugular vein, PU 2Fr 8cm, collar @ 1.1cm, round tip, VABM1B/25 attached, injector

Unit: each

Pouch Desc: Mouse JVC 2Fr PU 8cm VAB

Pouch: 3x7.5

Approx Dead Vol: 15 μ L

Connects to: 25ga

EO Sterilized: ☒

Intended Species: mouse

Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 8 cm

Material (IDxOD): PU 2Fr (.43x.69 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end

	Location	Fixed
Feature A: Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 1.1 cm	☐
Feature B:	b: cm	☐
Feature C:	c: cm	☐

Proximal end

Feature Z: Mouse VAB, 1ch magnetic, 25ga, bonded to catheter	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1: PNP3M injector

Part 2: