


Part Number: C20PU-MJV1825VAB



not to scale

Full Description: Catheter for mouse jugular vein, PU 2Fr 8cm, collar @ 1.25cm, round tip, VABM1B/25 attached, injector

Unit: each **Pouch Desc:** Mouse JVC 2Fr PU 8cm VAB **Pouch:** 3x7.5

Approx Dead Vol: 15 μ L **Connects to:** 25ga **EO Sterilized:** ☒

Intended Species: mouse **Vessel:** jugular vein

Usage Notes: collar position +/-0.5mm

TUBING

Tip: rounded **Length 1** **Length 2**

L1: 8 cm **L2:** cm

Material (IDxOD): PU 2Fr (.43x.69±.03mm) **Material (IDxOD):**

FEATURES

Distal (tip) end

Feature	Description	Location	Fixed
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 1.3 cm	☐
Feature B:		b: cm	☐
Feature C:		c: cm	☐

Proximal end

Feature	Description	Location	Fixed
Feature Z:	Mouse VAB, 1ch magnetic, 25ga, bonded to catheter	z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1: PNP3M injector **Part 2:**