


Part Number: C20PU-MJV2011

*not to scale*

Full Description: Mouse JVC, PU 2Fr 10cm, 11mm insertion length (a1.0cm), round tip. LS25 and SP25/12. Fits 25ga.

Unit: each Pouch Desc: Mouse JVC 2Fr PU L10cm 11mm Pouch: 3x7.5

Approx Dead Vol: 15 μ L Connects to: 25ga EO Sterilized: ☒

Intended Species: mouse Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 10 cm

Length 2

L2: cm

Material (IDxOD): PU 2Fr (.43x.69 $\pm$.03mm)

Material (IDxOD):

FEATURES

Distal (tip) end

Feature A: Small silicone collar for 2Fr (0.9mm OD), 1mm wide

Location

a: 1.0 cm

Fixed



Feature B:

b: cm



Feature C:

c: cm



Proximal end

Feature Z:

z: cm



Feature Y:

y: cm



Feature X:

x: cm

**PARTS ADDED TO POUCH**

Part 1: LS25

Part 2: SP25/12