


Part Number: C20PU-MJV2012



not to scale

Full Description: Mouse JVC, PU 2Fr 10cm, 12mm insertion length (a1.1cm), round tip. LS25 and SP25/12. Fits 25ga.

Unit: each

Pouch Desc: Mouse JVC 2Fr PU L10cm 12mm

Pouch: 3x7.5

Approx Dead Vol: 15 μL

Connects to: 25ga

EO Sterilized: ☒

Intended Species: mouse

Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 10 cm

Material (IDxOD): PU 2Fr (.43x.69±.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 1.1 cm	☒
Feature B:		b: cm	☐
Feature C:		c: cm	☐
Proximal end			
Feature Z:		z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1: LS25

Part 2: SP25/12