


Part Number: **C20PU-RFV1827**



not to scale

Full Description: Rat femoral vein catheter, 2Fr, 60cm, collar at 4cm, square tip, fits 25ga

Unit: each **Pouch Desc:** Rat FVC 2Fr PU 60cm **Pouch:** 4x10

Approx Dead Vol: 90 μ L **Connects to:** 25ga **EO Sterilized:**

Intended Species: rat **Vessel:** femoral vein

Usage Notes: non-survival; for use in MRI

TUBING

Tip: square **Length 1** L1: 60 cm **Length 2** L2: cm
Material (IDxOD): PU 2Fr (.43x.69 \pm .03mm) **Material (IDxOD):**

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 4.0 cm	<input type="checkbox"/>
Feature B:		b: cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**