


Part Number: C20PU-RSI2416



not to scale

Full Description: Catheter for rat small intestine, PU 2Fr 25cm, collars @ 0.7, 1.0cm. Fits 25ga.

Unit: each

Pouch Desc: Rat SIC 2Fr PU 25cm

Pouch: 4x10

Approx Dead Vol: 40 μL

Connects to: 25ga

EO Sterilized: ☒

Intended Species: rat

Vessel: small intestine

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 25 cm

Material (IDxOD): PU 2Fr (.43x.69±.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	a: 0.7 cm	☐
Feature B:	Small silicone collar for 2Fr (0.9mm OD), 1mm wide	b: 1.0 cm	☐
Feature C:		c: cm	☐

Proximal end			
Feature Z:		z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: