


Part Number: C30PU-BFV2208



not to scale

Full Description:

Catheter for rabbit femoral v/a, PU 3Fr 60cm, fixed collars @ 17, 17.5cm, rounded tip, sleeve. Fits 22ga.

Unit:

each

Pouch Desc:

Rabbit FVC 60cm

Pouch:

4x10

Approx Dead Vol:

190 μ L

Connects to:

22ga

EO Sterilized:

☒

Intended Species:

rabbit

Vessel:

femoral vein

Usage Notes:

TUBING

Tip:

rounded

Length 1

L1:

60

cm

Length 2

L2:

cm

Material (IDxOD):

PU 3Fr (.62x1.02 $\pm$.03mm)

Material (IDxOD):

FEATURES

		Location	Fixed
Distal (tip) end			
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 17.0 cm	☒
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 17.5 cm	☒
Feature C:		c: cm	☐
Proximal end			
Feature Z:	Silicone sleeve for 3Fr (1.7mm OD), 5mm wide	z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: