


Part Number: C30PU-BJV2408



not to scale

Full Description: Catheter for rabbit jugular vein, PU 3Fr 30cm, collars @ 5.0, 5.5cm, rounded tip. Fits 22ga.

Unit: each

Pouch Desc: Rabbit JVC 3Fr PU 30cm

Pouch: 4x10

Approx Dead Vol: 90 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: rabbit

Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 30 cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 5.0 cm	☐
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 5.5 cm	☐
Feature C:		c: cm	☐

Proximal end			
Feature Z:		z: cm	☐
Feature Y:		y: cm	☐
Feature X:		x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: