


Part Number: C30PU-GFV2421



not to scale

Full Description: Catheter for guinea pig femoral vein, PU 3Fr 30cm, collars @ 7.0, 7.5cm, rounded tip. Fits 22ga.

Unit: each

Pouch Desc: GP FVC 3Fr PU 30cm

Pouch: 4x10

Approx Dead Vol: 90 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: guinea pig

Vessel: femoral vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 30 cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end

	Location	Fixed
Feature A: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 7.0 cm	☐
Feature B: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 7.5 cm	☐
Feature C:	c: cm	☐

Proximal end

Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: