


Part Number: C30PU-GJV2306



not to scale

Full Description: Catheter for GP jugular vein, PU 3Fr 14cm, collar @ 2.0cm, rounded tip. Fits 22ga.

Unit: each

Pouch Desc: Rat JVC 3Fr PU 14cm a2.0

Pouch: 3x7.5

Approx Dead Vol: 40 μL

Connects to: 22ga

EO Sterilized: ☒

Intended Species: guinea pig

Vessel: jugular vein

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 14 cm

Material (IDxOD): PU 3Fr (.62x1.02±.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end

	Location	Fixed
Feature A: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 2.0 cm	☐
Feature B:	b: cm	☐
Feature C:	c: cm	☐

Proximal end

Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: