


Part Number: **C30PU-PJV1805**



not to scale

Full Description: Catheter for minipig jugular vein or carotid artery, PU 3Fr 30cm, collar @ 9cm. Fits 22ga.

Unit: each **Pouch Desc:** Minipig JVC 3Fr PU 30cm **Pouch:** 4x10

Approx Dead Vol: 95 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: minipig **Vessel:** jugular vein

Usage Notes:

TUBING

Tip: rounded **Length 1** L1: 30 cm **Length 2** L2: cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 9.0 cm	<input type="checkbox"/>
Feature B:		b: cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**