


Part Number: **C30PU-RBD1424**



not to scale

Full Description: Catheter for rat bile duct/duodenum, PU 3Fr 40cm, round and beveled tips

Unit: each **Pouch Desc:** Rat BDC 3Fr PU 40cm **Pouch:** 4x10

Approx Dead Vol: 130 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** bile duct/duodenum

Usage Notes: Round tip end intended for duodenum, beveled end for bile duct. Cut in middle for sampling. Can be attached to Instech bile collection harness (VAHD115AB, VAHD115L).

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 40 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 0.3 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 0.6 cm	<input type="checkbox"/>
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 5.0 cm	<input type="checkbox"/>
Proximal end			
Feature Z:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	z: 1.0 cm	<input type="checkbox"/>
Feature Y:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	y: 5.0 cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**