


Part Number: **C30PU-RBD1801**



*not to scale*

**Full Description:** Catheter for rat bile duct/duodenum, PU 3Fr 40cm, round and beveled tips

**Unit:** each **Pouch Desc:** Rat BDC 3Fr PU 40cm **Pouch:** 4x10

**Approx Dead Vol:** 130  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** bile duct/duodenum

**Usage Notes:** Round tip end intended for duodenum, beveled end for bile duct. Cut in middle for sampling. Can be attached to Instech bile collection harness (VAHD115AB, VAHD115L).

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 40 cm **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm) **Material (IDxOD):** cm

## FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 0.3 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 0.6 cm	<input type="checkbox"/>
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 5.0 cm	<input type="checkbox"/>
Proximal end			
Feature Z:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	z: 1.0 cm	<input type="checkbox"/>
Feature Y:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	y: 5.0 cm	<input type="checkbox"/>
Feature X:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	x: 8.0 cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**