


Part Number: C30PU-RBL2401



not to scale

Full Description: Catheter for rat bladder, PU 3Fr 60cm, collars @ 0.2, 0.5, 4.0cm, rounded tip. Fits 22ga.

Unit: each

Pouch Desc: Rat bladder cath 3Fr 60cm

Pouch: 4x10

Approx Dead Vol: 190 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: rat

Vessel: bladder

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 60 cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end

	Location	Fixed
Feature A: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 0.2 cm	☐
Feature B: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 0.5 cm	☐
Feature C: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 4.0 cm	☐

Proximal end

Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: