


Part Number: C30PU-RBL2407



not to scale

Full Description: Catheter for rat bladder, PU 3Fr 30cm, holes at 0.2, 0.4, collars @ 0.5, 0.8, 4.0cm, rounded tip

Unit: each

Pouch Desc: Rat bladder cath 3Fr 30cm

Pouch: 4x10

Approx Dead Vol: 90 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: rat

Vessel: bladder

Usage Notes:

TUBING

Tip: rounded

Length 1

L1: 30 cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Perfusion holes, 2 sets, 2mm apart	a: 0.2 cm	☐
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 0.5 cm	☐
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 0.8 cm	☐
Proximal end			
Feature Z:		z: cm	☐
Feature Y:		y: cm	☐
Feature X:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	x: 26.0 cm	☐

PARTS ADDED TO POUCH

Part 1:

Part 2: