


Part Number: **C30PU-RFV1303VAB**



*not to scale*

**Full Description:** Catheter for rat fem vein, PU 3Fr 19cm, collars @ 4.5, 7 and 16cm, VAB95BS attached

**Unit:** each      **Pouch Desc:** VAB w/ FVC 3Fr PU 19cm      **Pouch:** 4x10

**Approx Dead Vol:** 75  $\mu$ L      **Connects to:** 22ga      **EO Sterilized:**

**Intended Species:** rat      **Vessel:** femoral vein

**Usage Notes:**

## TUBING

**Tip:** rounded      *Length 1*      *Length 2*  
**L1:** 19 cm      **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm)      **Material (IDxOD):** cm

## FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
<b>Feature A:</b>	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	a: 4.5 cm	<input type="checkbox"/>
<b>Feature B:</b>	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	b: 7.0 cm	<input type="checkbox"/>
<b>Feature C:</b>		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
<b>Feature Z:</b>	Rat VAB, 22ga connector	z: cm	<input type="checkbox"/>
<b>Feature Y:</b>	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	y: 3.0 cm	<input type="checkbox"/>
<b>Feature X:</b>		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:** VAH6M injector      **Part 2:**