


Part Number: **C30PU-RFV1308**



not to scale

Full Description: Catheter for rat femoral vein, PU 3Fr 25cm, collars @ 3.0, 3.5cm. Fits 22ga.

Unit: each **Pouch Desc:** Rat FVC 3Fr PU 25cm **Pouch:** 4x10

Approx Dead Vol: 75 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** femoral vein

Usage Notes:

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 25 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):**

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 3.0 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 3.5 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**