


Part Number: **C30PU-RFV1408**



*not to scale*

**Full Description:** Catheter for rat femoral vein, PU 3Fr 19cm, fixed collars @ 3.7, 4.9 and 7.4cm. Fits 22ga.

**Unit:** each **Pouch Desc:** Rat FVC 3Fr PU 19cm **Pouch:** 4x10

**Approx Dead Vol:** 60  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** femoral vein

**Usage Notes:** collars in same position as PU-C35-BAY2

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 19 cm **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm) **Material (IDxOD):** cm

## FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	a: 3.7 cm	<input checked="" type="checkbox"/>
Feature B:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	b: 4.9 cm	<input checked="" type="checkbox"/>
Feature C:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	c: 7.4 cm	<input checked="" type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**