


Part Number: **C30PU-RFV1456**



not to scale

Full Description: Catheter for rat femoral vein, PU 3Fr 76cm, collars @ 2.5, 3.0cm. 15in VAHBPU-T22 with LS22.

Unit: each **Pouch Desc:** Rat FVC 3Fr PU 76cm kit **Pouch:** 4x10

Approx Dead Vol: 245 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** femoral vein

Usage Notes:

TUBING

Tip: rounded **Length 1** L1: 76 cm **Length 2** L2: cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):**

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 2.5 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 3.0 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: 15in VAHBPU-T22 with LS22 **Part 2:**