


Part Number: **C30PU-RFV1604**



*not to scale*

**Full Description:** Catheter for rat femoral vein, PU 3Fr 30cm, collars @ 4.5, 5.0cm. Fits 22ga.

**Unit:** each **Pouch Desc:** Rat FVC 3Fr PU 30cm **Pouch:** 4x10

**Approx Dead Vol:** 90  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** femoral vein

**Usage Notes:**

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 30 cm **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm) **Material (IDxOD):**

## FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
<b>Feature A:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 4.5 cm	<input type="checkbox"/>
<b>Feature B:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 5.0 cm	<input type="checkbox"/>
<b>Feature C:</b>		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
<b>Feature Z:</b>		z: cm	<input type="checkbox"/>
<b>Feature Y:</b>		y: cm	<input type="checkbox"/>
<b>Feature X:</b>		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**