


Part Number: **C30PU-RFV1802**



not to scale

Full Description: Catheter for rat femoral vein, PU 3Fr 76cm, collars @ 3.7, 4.9, 7.4cm. Fits 22ga.

Unit: each **Pouch Desc:** Rat FVC PU 3Fr 76cm **Pouch:** 4x10

Approx Dead Vol: 250 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** femoral vein

Usage Notes:

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 76 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 3.7 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 4.9 cm	<input type="checkbox"/>
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 7.4 cm	<input type="checkbox"/>

<i>Proximal end</i>		<i>Location</i>	<i>Fixed</i>
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**