


Part Number: **C30PU-RFV1850VAB**



not to scale

Full Description: Catheter for rat fem vein, PU 3Fr 22cm, collars @ 5.0, 9.0cm, VABR1B/22 attached

Unit: each **Pouch Desc:** VAB FVC 3Fr L22 a5.0cm **Pouch:** 4x10

Approx Dead Vol: 75 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** femoral vein

Usage Notes: bead locations +/-1mm, length +/-3mm; volume estimate +/-15%

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 22 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02±.03mm) **Material (IDxOD):** cm

FEATURES

Distal (tip) end

Feature	Description	Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 5.0 cm	<input checked="" type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 9.0 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>

Proximal end

Feature Z:	Rat VAB, 1ch, 22ga (VABR1B/22)	z: cm	<input checked="" type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: PNP3M injector **Part 2:**