


Part Number: **C30PU-RFV1858VAB**



not to scale

**Full Description:** Catheter for rat fem vein, PU 3Fr 22cm, collars @ 5.8, 9.8cm, VABR1B/22 attached

**Unit:** each **Pouch Desc:** VAB FVC 3Fr L22 a5.8cm **Pouch:** 4x10

**Approx Dead Vol:** 75  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** femoral vein

**Usage Notes:** bead locations +/-1mm, length +/-3mm; volume estimate +/-15%

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 22 cm **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02±.03mm) **Material (IDxOD):** cm

## FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 5.8 cm	<input checked="" type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 9.8 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:	Rat VAB, 1ch, 22ga (VABR1B/22)	z: cm	<input checked="" type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:** PNP3M injector **Part 2:**