


Part Number: **C30PU-RFV1863VAB**



not to scale

Full Description: Catheter for rat fem vein, PU 3Fr 22cm, collars @ 6.3, 10.3cm, VABR1B/22 attached

Unit: each **Pouch Desc:** VAB FVC 3Fr L22 a6.3cm **Pouch:** 4x10

Approx Dead Vol: 75 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** femoral vein

Usage Notes: bead locations +/-1mm, length +/-3mm; volume estimate +/-15%

TUBING

Tip: rounded *Length 1* *Length 2*
 L1: 22 cm L2: cm
Material (IDxOD): PU 3Fr (.62x1.02±.03mm) **Material (IDxOD):** cm

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 6.3 cm	<input checked="" type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 10.3 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
Feature Z:	Rat VAB, 1ch, 22ga (VABR1B/22)	z: cm	<input checked="" type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: PNP3M injector **Part 2:** cm